

MOBILE VACCINE REFRIGERATED CART



COVID-19 VACCINATION REFRIGERATED CART

VR4300K-AC

Features:

- Allows for extended storage for vaccines prior to dispensing
- Refrigerator keeps vaccines down to 2° C
- Lab and medical-grade quality 2.4 cubic feet of refrigerated capacity
- Interior refrigerator dims are 20.5" H x 15.5" W x 13.5" D
- User-programmable, audible and visual high/low alarms
- Hospital grade power strip and 15 foot long cord included
- Secured drawer for storage of syringes masks, gauze and other vaccination supplies
- Key lock on drawers and refrigerator
- Removable plastic top
- Fully welded heat sealed powder coat body and drawers with 12 year warranty

Colors:

Shown in Hammertone Blue, but available in more than twenty color options. View them online: www.harloff.com/finish-options/

Shown with Refrigerator, also available in Freezer models VF4300K-AC & VF4250K-AC



Mobile Vaccine Refrigerated Cart for the Application of COVID-19 Vaccinations, Three Drawer, Key Lock, No Accessories

Engineered & Manufactured in the U.S.



Warranty:

12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Overall Cart Dimensions:

- Internal Drawer Dimensions: 23" W x 17" D (58.42 cm W x 43.18 cm D)
- Cart dims: 51.5" H x 37.5" W x 26" D (130.8 cm H x 95.3cm W x 66cm D)
- Cart net weight: 145 lbs (65.8 kg)

Drawer Configuration:

- (3) 3" (7.6 cm) and (1) Compartment with pre-installed refrigerator

Optional Accessories:



Glove Box Holder Triple Direct Mount



Chart/document holders



Sharps container & mounting bracket



Waste container & mounting bracket



Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300