



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Wood Vinyl 294 Capacity Unit Dose Medication Cart, Standard Package #WV294BOX-CM

This vinyl coated aluminum punch card cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of Long Term Care use.



Skirt Option



Cherry
 Mahogany
 (CM)

WV294BOX-CM

Color Options:

Shown in Cherry Mahogany.

Harloff offers two standard vinyl color options, plus additional special order vinyl color options with a minimum cabinet purchase and extended lead time.

Specifications

Construction:

- Vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate

Features:

- Light weight vinyl coated aluminum
- Single column medication cart with 294 unit dose capacity
- Key locking cabinet with locking narcotics box located in bottom drawer
- Adjustable dividers for unit dose boxes and liquid storage
- Ball-bearing, full extension drawer slides
- Pull-out shelf
- 5" ball-bearing premium rubber casters
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture

Skirt Option:

- #WV-3WSKIRT-CM – Skirt option to cover wheels

Drawer Configuration:

- 1 ea – 3" multipurpose storage drawer
- 6 ea – 3" multipurpose storage drawers with adjustable dividers
- 1 ea - 9.75" punch card drawer with adjustable divider and narcotics box

Assembled Cart Dimensions:

- 45" H x 34.5" W x 24.5" D
- Writing Surface: 44" (111.8)
- Cart weight: 145 lbs (65.9 kg)
- Approx shipping weight: 205 lbs (93.2 kg)
- Cart will ship on a pallet