



*Wood Vinyl*

**WV3256K-CM**

## TALL MINI WIDE ANESTHESIA CART

Wood Vinyl Covered Aluminum with a Key Lock

### Features:

- Lightweight vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate
- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Key lock with two keys
- Bumper system extends wheelbase and stabilizes cart
- Ball-bearing, full extension drawer slides
- 5" ball-bearing full swivel casters – one brake, one directional
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture
- Latex-free
- Twelve Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



**Engineered & Manufactured in the U.S.**



(CM) Cherry Mahogany

### Skirt Option:

- #WV-2WSKIRT-CM – Skirt option to cover wheels

### Colors:

Shown in Cherry Mahogany. Harloff offers two standard vinyl color options, plus additional special order vinyl color options with a minimum cabinet purchase and extended lead time.



#### Drawer Configuration:

- 12.6" W x 16.75" D (32cm x 42.5cm)
- Three 3.1" drawers (7.9 cm), two 6.4" drawers (16.4 cm) and one 9.8" drawer (25 cm)

#### Overall Cart Dimensions:

- 43.75" H x 19.5" W x 25" D (110.5 cm H x 49.5 cm W x 63.5 cm D)
- Cart weight: 97 pounds (44 kg)
- Approx. shipping weight: 107 pounds (48.5 kg)

Harloff manufactures a complete line of medical furniture for use in Acute Care Facilities. These products provide an upscale look when compared to typical metal options. Our product line includes isolation, emergency, treatment and medication carts.



View similar products online at

<https://www.harloff.com/product-category/medical-furniture-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300