



Wood Vinyl

WV540PC-DW

PUNCH CARD MEDICATION CART

540 Capacity Wood Vinyl Covered Aluminum

Features:

- Lightweight vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate
- Single column medication cart with 540 punch card capacity
- Key locking cabinet and narcotics box located in top drawer
- Adjustable row dividers for punch card and liquid storage
- Ball-bearing, full extension drawer slides
- Pull-out shelf
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture
- Twelve Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Drawer Configuration:

- 1 ea – 3" multipurpose storage drawer
- 1 ea – 9.75" punch card drawer with locking narcotics box and adjustable dividers
- 2 ea - 9.75" punch card drawer with adjustable dividers

Skirt Option:

- #WV-3WSKIRT-DW – Skirt option to cover wheels

Colors:

Shown in Driftwood Gray color. Harloff offers standard vinyl color options, plus additional special order vinyl color options with a minimum cabinet purchase and extended lead time.



(DW)
Driftwood Gray

This vinyl coated aluminum punch card cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of Long Term Care use.

Overall Cart Dimensions:

- 45" H x 34.5" W x 24.5" D
- Writing Surface: 44" (111.8)
- Cart weight: 127 lbs (57.7 kg)
- Approximate shipping weight: 157 lbs (71.2 kg)



Engineered & Manufactured in the U.S.



View similar products online at

<https://www.harloff.com/product-category/medication-furniture-carts-cabinets/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com

www.harloff.com • Ph: (719) 637-0300