



# Wood Vinyl

## WV6200-CM

### TALL ALUMINUM TREATMENT CART Wood Vinyl with Storage Compartment & Key Lock

#### Features:

- Lightweight vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate
- 21" Storage compartment, locking doors and adjustable shelf
- Key lock on 2 top drawers with 2 keys, additional key lock on storage compartment with 2 keys
- Ball-bearing, full extension drawer slides
- Cart is washable with soap and water
- 5" (12.7cm) ball-bearing premium casters, all full swivel, one braking, one directional
- This cart is designed to look like a piece of furniture
- Twelve Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



Engineered & Manufactured in the U.S.

#### Skirt Option:

- #WV-3WSKIRT-CM – Skirt option to cover wheels

#### Colors:

Shown in Cherry Mahogany. Harloff offers two standard vinyl color options, plus additional special order vinyl color options with a minimum cabinet purchase and extended lead time.



(CM) Cherry Mahogany



#### Drawer Configuration:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)
- One 3.1" drawer (7.9 cm), one 6.4" drawer (16.4 cm) and one 21" storage compartment with adjustable internal shelf

#### Overall Cart Dimensions:

- 41.5" H x 32" W x 22" D (105.4 cm H x 81.3 cm W x 55.9 cm D)
- Cart weight: 100 pounds (45.4 kg)
- Shipping dims: 48" H x 28" W x 39" D (121.9 cm H x 71.1 cm W x 99.1 cm D)
- Approx. shipping weight: 120 pounds (54.4 kg)

Bring the warm appearance of wood to your treatment and procedure carts with vinyl-coated aluminum carts. Designed and built to support culture change initiatives, this cart gives the buyer the best of both worlds, combining the aesthetics of wood with the long-term functionality of a Harloff metal cart.



View similar products online at <https://www.harloff.com/product-category/medical-furniture-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300