



Wood Vinyl

WV6350-CM

SHORT HEIGHT TREATMENT CART

Wood Vinyl Covered Aluminum with a Key Lock

Features:

- Lightweight vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate
- Pull out shelf
- Key lock with 2 keys
- Ball-bearing, full extension drawer slides
- 5" (12.7cm) ball-bearing premium casters, all full swivel, one braking, one directional
- This cart is designed to look like a piece of furniture
- Cart is washable with soap and water
- Twelve Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



Engineered & Manufactured in the U.S.

Skirt Option:

- #WV-3WSKIRT-CM – Skirt option to cover wheels

Colors:

Shown in Cherry Mahogany.

Harloff offers two standard vinyl color options, plus additional special order vinyl color options with a minimum cabinet purchase and extended lead time.



(CM) Cherry Mahogany



Drawer Configuration:

- Two 3.1" drawers (7.9 cm), one 6.4" drawer (16.4 cm) and one 13.1" drawer (33.3 cm)

Overall Cart Dimensions:

- 38.5" H x 30" W x 22.5" D (97.8 cm H x 76.2 cm W x 57.2 cm D)
- Writing Surface: 38" height (96.5 cm)
- Cart weight: 110 pounds (49.9 kg)
- Approx. shipping weight: 140 pounds (63.5 kg)

Bring the warm appearance of wood to your treatment and procedure carts with vinyl-coated aluminum carts. Designed and built to support culture change initiatives, this cart gives the buyer the best of both worlds, combining the aesthetics of wood with the long-term functionality of a Harloff metal cart.

View similar products online at

<https://www.harloff.com/product-category/medical-furniture-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com

www.harloff.com • Ph: (719) 637-0300

HARLOFF
INNOVATING STORAGE FOR CARE