

Toll: (800) 433-4064 Ph: (719) 637-0300 medicalsales@harloff.com www.harloff.com

MR-Conditional Line/ Model: MR6K

Aluminum Six Drawer Anesthesia Cart, Key Lock, Standard Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

 Non-magnetic or "weekly magnetic" materials including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware

Features:

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; two swivel
 two with brake, one directional
- · Painted aluminum breakaway locking system
- BEST® brand lock on cabinet, with 2 keys
- Protective corner bumpers
- · Painted aluminum pull-out shelf
- Antimicrobial plastic top with integrated push handles, removable
- Cart is delivered fully assembled with pre-threaded accessory mounting holes
- MR-Conditional accessories available
- Latex free
- · Durable, hard-baked powder coat finish

Drawer Configuration:

- (1) 4" (10.2 cm)
- (2) 3.25" (8.3 cm)
- (2) 6.5" (16.5 cm)
- (1) 9.75" (24.8 cm)

Internal Drawer Configuration:

• 23.9" W x 17" D (60.8 cm x 43.2 cm)

Overall Cart Dimensions:

- 45.75" H x 30.5" W x 23.75" D (116.6 cm H x 77.5 cm W x 60.3 cm D)
- Writing surface: 43.5" (110.5cm)
- Cart net weight: 108 lbs (49.5 kg)
- Approximate shipping weight: 146 lbs (66.4 kg)
- Approximate shipping dims: 55" H x 37" W x 26" D (139.7cm H x 94cm W x 66cm D)