

# THE CHOICE IS YOURS WITH A-SERIES AND M-SERIES CARTS!

CHOOSE YOUR HEIGHT, WIDTH AND DRAWER CONFIGURATION  
CHOOSE YOUR LOCK - YOUR COLOR - YOUR ACCESSORIES

## STEEL OR ALUMINUM?

*The A-Series is our line of ultra light aluminum carts with optimal curb appeal. Enjoy top-quality innovative storage and organization with a sleek, modern appearance. The A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function.*

*The M-Series is our line of durable steel carts designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top help make things easier for your facility.*



Engineered and  
Manufactured in the U.S.

### A-SERIES FEATURES:

- Lightweight aluminum fully welded uni-body construction
- Curved drawer pulls, durable and elegant
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Ergonomic side push handles can be mounted on either side
- Antimicrobial plastic top, easily removable and replaceable
- Fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered fully assembled



MPA3030E06

### M-SERIES FEATURES:


- Choice of lock options
- Full extension drawers slides for easy access
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Durable scratch resistant hard baked powder coat finish
- Ergonomic side push handles can be mounted on either side
- Industry best 5" quiet rolling Tente® casters - 1 steering caster & 2 braking casters
- Body has pre-threaded accessory mounting holes
- Antimicrobial, removable plastic top



MDS3030B16

# M-SERIES/A-SERIES CARTS


## THREE WIDTHS TO CHOOSE FROM:




## CHOOSE YOUR DRAWER CONFIGURATION:

1. Choose a cart size (X-Short, Short, Medium and Tall)
2. Select any combination of 3", 6", 9" or 12" drawers that fit within your cart's configurable drawer space.


5 drawers x 6" height = 30"




3 drawers x 9" height = 27"



2 drawers x 12" height = 24"



7 drawers x 3" height = 21"



## CHOOSE YOUR LOCK:

- B - Breakaway Lock • K - Key-Lock • E - Basic Electronic Pushbutton Lock • EKC - Electronic Lock w/Keypad
  - ELP - Electronic Lock w/Keypad and Proximity Reader • WKC - Wireless Electronic Lock w/Keypad • WELP - Wireless Electronic Lock w/Keypad and Proximity Reader
- View locking options online at: <https://www.harloff.com/locking-options/>

## CHOOSE YOUR COLOR(S):

Available in over 20 color options. An anti-microbial finish is available in Lemon Drop Yellow. View all color options online at: <https://www.harloff.com/finish-options/>

## CHOOSE YOUR ACCESSORIES:

A full line of drawer dividers and exterior accessories are available for all cart sizes. View accessory options online at: <https://www.harloff.com/accessories/>

## COMMON DRAWER CONFIGURATIONS:



MDS3030K04

1-3, 3-9

MDS3030K05

2-3, 1-6  
2-9

MDS3030K06

3-3, 2-6  
1-9

MDS3030K16

4-3, 1-6  
1-12



MDS3024K04

1-3, 2-6  
1-9

MDS3024K05

3-3, 1-6  
1-9

MDS3024K15

2-3, 3-6



MDS3027K04

1-3, 1-6  
2-9

MDS3027K05

2-3, 2-6  
1-9

MDS3027K06

4-3, 1-6  
1-9

MDS3027K16

3-3, 3-6



MDS3021K04

1-3, 3-6

MDS3021K14

2-3, 1-6  
1-9

MDS3021K05

3-3, 2-6

# MPA30 CARTS



## SHORT A-SERIES

**Configurable Drawer Space: 24"**

**37.25" H x 36.75" W x 22" D**

36.125" Height of  
Work Surface

**X-SHORT A-SERIES**  
**Configurable Drawer Space: 21"**

**34" H x 36.75" W x 22" D**

32.875" Height of  
Work Surface



MPA3021K04



MPA3024E05

**MEDIUM A-SERIES**  
**Configurable Drawer Space: 27"**

**40.5" H x 36.75" W x 22" D**

39.375" Height of  
Work Surface



MPA3027B06

**TALL A-SERIES**  
**Configurable Drawer Space: 30"**  
**43.75" H x 36.75" W x 22" D**  
42.625" Height of  
Work Surface



MPA3030EKC06

# MDS30 CARTS



## SHORT M-SERIES

**Configurable Drawer Space: 24"**

**37.25" H x 36.75" W x 22" D**

36.125" Height of  
Work Surface

**X-SHORT M-SERIES**  
**Configurable Drawer Space: 21"**

**34" H x 36.75" W x 22" D**

32.875" Height of  
Work Surface



MDS3021E05



MDS3024E05

**MEDIUM M-SERIES**  
**Configurable Drawer Space: 27"**

**40.5" H x 36.75" W x 22" D**

39.375" Height of  
Work Surface



MDS3027ELP06

**TALL M-SERIES**  
**Configurable Drawer Space: 30"**  
**43.75" H x 36.75" W x 22" D**  
42.625" Height of  
Work Surface



MDS3030B06



# MPA24 AND MPA18 CARTS

## MPA2430EKC16

TALL

30" of Configurable  
Drawer Space  
43.75" H x 30" W x  
22" D, 42.625" Height  
of Work Surface



## MPA2427E05

MEDIUM

27" of Configurable  
Drawer Space  
40.5" H x 30" W x 22" D  
39.375" Height of  
Work Surface



## MPA2421K14

X-SHORT

Configurable Drawer  
Space: 21",  
34" H x 30" W x 22" D  
32.875" Height of  
Work Surface



SHORT

Configurable Drawer  
Space: 24", 37.25" H x 30" W x 22" D  
36.125" Height of Work Surface

## MPA2424B05

## M3PA1830K06

TALL

30" of Configurable  
Drawer Space  
43.75" H x 24.75" W x 22" D  
42.625" Height of  
Work Surface



## MPA1827E05

MEDIUM

27" of Configurable  
Drawer Space  
40.5" H x 24.75" W x 22" D  
39.375" Height of  
Work Surface



## MPA1821EKC14

X-SHORT

Configurable  
Drawer  
Space: 21",  
34" H x 24.75"  
W x 22" D  
32.875"  
Height of  
Work Surface



SHORT

Configurable  
Drawer Space: 24",  
37.25" H x 24.75" W x 22" D  
36.125" Height of Work Surface

## MPA1824L03

## MDS2421K14

### X-SHORT

Configurable

Drawer Space: 21",  
34" H x 30" W x 22" D

32.875" Height of  
Work Surface



### SHORT

Configurable

Drawer Space: 24",  
37.25" H x 30" W x 22" D

36.125" Height of Work Surface

## MDS2424E05

# MDS24 AND MDS18 CARTS

## MDS2427EKC05

### MEDIUM

27" of Configurable  
Drawer Space  
40.5" H x 30" W x 22" D

39.375" Height of  
Work Surface



## MDS2430KC5-18EH

### TALL

30" of Configurable  
Drawer Space  
43.75" H x 30" W x 22" D  
42.625" Height of Work Surface



Cart is  
shown with an  
equipment holder  
upgrade!

## MDS1821E14

### X-SHORT

Configurable

Drawer  
Space: 21",  
34" H x 24.75"  
W x 22" D

32.875" Height  
of Work Surface



### SHORT

Configurable

Drawer Space: 24",  
37.25" H x 24.75" W x 22" D

36.125" Height of Work Surface

## M3DS1824K03

## MDS1827WKC05

### MEDIUM

27" of Configurable  
Drawer Space  
40.5" H x 24.75" W x  
22" D, 39.375" Height  
of Work Surface



## MDS1830EKC06

### TALL

30" of Configurable  
Drawer Space  
43.75" H x 24.75" W x 22" D  
42.625" Height of  
Work Surface



# M-SERIES/A-SERIES OPTIONS

## EQUIPMENT HOLDERS

Replace 9", 12" or 18" of vertical drawer space at the top or bottom of your M-Series or A-Series cart with an open or closed door equipment holder, perfect for storing large medical equipment. Equipment holders feature a 17" deep holding area that replaces standard drawer areas and comes with a 2" hole in the rear for cords and cables to pass through.

**MDS3030KC4+MD30-DRW9-EH**

**MDS2430KC4+MD24-DRW18-EHDR**



**MDS3030B16+MD30-FDR12**



## DRAWER OPTIONS

For carts with storage areas, you have several options available to you for these enclosures including flip up doors, locking double doors and adjustable shelves. Don't see what you want? Harloff can create custom drawer designs to match unique requirements.

## BUMPERS AND CASTERS

Standard M-Series and A-Series carts come equipped with fully curved "pontoon-style" protective bumper guards and premium 5" quiet rolling Tente® casters, but other options are available if you have (for example) height or width space limitations or other considerations. We offer two non-custom bumper options (smaller gray bumpers or no bumpers) and two caster options (3" plate casters or 3" stem casters).

**MDS1824E04G3**

(gray bumpers, 3" plate casters)

**MDS2421E05N5**

(no bumpers, 5" plate casters)

**M3DS1830K05**

(no bumpers, 3" stem casters)



## CUSTOM OPTIONS

M-Series/A-Series carts are highly adaptable, customizable and can be used for a variety of medical purposes. Let us help build the perfect custom cart for you!

**MPA3030EC5**

Includes an open tray holder option MD30-DRW15-EH-HH





Easily add the look of wood to any M-Series or A-Series cart by adding vinyl-coated aluminum drawer fronts. This combines the aesthetics of furniture with the ability to customize and the long-term functionality of a Harloff metal cart.

## WOOD VINYL DRAWER FRONTS

### MDS1830KC7-CM

This mini-width M-Series cart features seven drawers with Cherry Mahogany wood vinyl drawer fronts.

### MDS3030EC6-AN

Give your M-Series cart a home-like feel with American Natural drawer fronts.

### MPA3024EC5-AN



Wood vinyl color options include American Natural, Cherry Mahogany and several other types of wood grain colors.



## WOOD GRAIN FINISH CARTS

### WV3256K-CM

We offer full wood grain finish options for several cart styles including crash carts (WV6400B) anesthesia carts (WV3256K), and general treatment carts (WV6200, WV6350, WV6450, WV6520, WV6400SC). Each

"medical furniture" cart features the same full extension drawer slides and 5" premium ball-bearing casters of our standard carts, but comes with a durable vinyl coated aluminum body, a pull out shelf and gray bumpers. These carts are ideal for facilities supporting Culture Change initiatives by giving carts a home-like feel.

### WV6200-CM

Similar to the MDS3030K02-21DR, this wood vinyl coated aluminum cart has a 21" tall bottom storage area.



### WV6400B-CM



# 10 REASONS TO CHOOSE HARLOFF ACUTE CARE CARTS

## 1 TENTE, World Leader in Casters

Tente is renowned for having the **Highest Quality Casters** in the industry, providing the best longevity and durability.



## 6 Proudly Made in the USA

Harloff products are proudly **Designed, Engineered and Manufactured in the United States** at our Colorado Springs, CO facility.



## 2 Antimicrobial Silver Infused Plastic Tops

This high touch area employs silver technology to enhance **Infection Control**. The replaceable top is easily removed for thorough cleaning.



## 7 Pair of Ergonomic Push Handles Included

Ergonomic **side handles** are **included** on every M-Series and A-Series model (excluding narrow carts) and can be mounted on either side of the cart.

## 3 Radius Corners and Flat Sides

Smooth surfaces assist **Infection Control** compliance by providing an easy surface to clean and sanitize. Difficult to clean corner channels are not used on Harloff carts.



## 8 Drawer Customizations/Options

**"One Size Does Not Fit All!"** Harloff will design and build your cart to meet your needs. There are **No Additional Charges** for changes to drawer configuration and accessory placement.

## 4 Heat Sealed Powder Coat Painted Inside and Out

We use a multi-step process that provides a hard, scratch resistant surface for extreme durability. **Harloff's Non Porous Surface** is easy to sanitize and enhances infection control protocols.



## 9 Pontoon Bumpers

Our **Large Bumpers** **Provide Protection** for the back, front and sides of the cart and protect facility walls from damage. Other bumpers are also available.

## 5 12 Year Warranty

Our **Manufacturer's Warranty** is **Unmatched** in the industry. We build our carts to last and stand behind our work.



## 10 Color Options

With a selection of over **20 Colors** to choose from, you can always get just the color you need. You can even mix and match the cart body and drawers with two separate colors.





# QUICK-RESPONSE CARTS DESIGNED FOR EMERGENCY ENVIRONMENTS

Harloff provides emergency carts to meet your unique needs. In addition to the many standard configurations shown, we are well known for our custom capabilities. Whether it's changing carts to custom locks, replacing drawer configurations, creating specialized accessories for unique facility requirements or designing and building new specially-sized carts, just ask our Sales Team for a quote to meet your specific needs.

MDS3030B06+MD30-EMG3

## EMERGENCY CRASH CARTS

Emergency carts, code carts and crash carts are designed for storing supplies needed in emergency situations. Harloff crash carts feature quality construction, multiple configurations and a 12 year warranty on manufactured parts.

Harloff is the market leader in storage solutions designed for emergency departments, ICU coronary care and urgent care centers. Our crash carts are used in world-class hospitals and medical facilities around the world.



Engineered and  
Manufactured in the U.S.



Includes  
the MD30-EMG3  
comprehensive Super  
Stat accessory  
package!

# EMERGENCY CRASH CARTS

## MDS3030B06+MD30-EMG3

Replaces  
the 6411

This top-of-the-line Harloff emergency cart has plenty of drawer space for supplies, the mobility to arrive quickly at any emergency and the accessories required to hold necessary medical equipment. Includes Harloff's MDS30-EMG3 Super Stat Emergency Package, a comprehensive array of crash cart accessories to keep your facility prepared for any contingency, plus a breakaway bar with 100 seals.



Equip your cart with an accessory package for the ideal emergency crash cart setup, or select from other options! M-Series and A-Series offer you choices.

## MD30-EMG3 ACCESSORY PACKAGE

For MDS30/MPA30 carts, this Super Stat Package includes:

- Cardiac board and brackets (MD-CARDBRD)
- Oxygen tank holder (O2HLDR)
- IV pole (MD-IV2)
- (2) Utility hooks (UHOOKS2)
- Large Utility Basket (BASKETDM)
- Chart Holder (CHARTHLDRDM)
- Toprail, matching cart color (MD30-TOPRAIL)
- Defibrillator Shelf (MD-DEFIB)
- Locking Sharps Container (5qt) and Glove Box Holder (LKSHRPHOLD5Q-GLV)
- Suction shelf (SUCTIONSHLF)
- (1) Standard 3" Drawer Dividers (MD30-DIV3-S)
- (1) Standard 6" Drawer Dividers (MD30-DIV6-S)
- (1) Standard 3" Drawer Tray Dividers (MD30-TRAYDIV3-S)



## MD30-EMG1 ACCESSORY PACKAGE

For MDS30/MPA30 carts, this Standard Package includes:

- IV pole (MD-IV2)
- Cardiac board and brackets (MD-CARDBRD)
- Oxygen tank holder (O2HLDR)
- Drop shelf (MD-DRPSHLF)
- (2) Utility hooks (UHOOKS2)
- (2 Sets) Standard 3" Drawer Dividers (MD30-DIV3-S)



## MDS3030B06+MD30-EMG1

Replaces  
the 6401

This M-Series tall emergency crash cart includes Harloff's MD30-EMG1 Emergency Package which features the most common requirements for medical crash carts. This includes an IV pole, oxygen tank holder, cardiac board, drop shelf, utility hooks and drawer dividers. Other features include industry best Tente® casters for quick movement, top quality drawer slides which provide supply access and a breakaway locking bar with 100 seals to ensure each cart is stocked and ready. Harloff crash carts are known and depended upon by thousands of medical professionals.



QuickShip

available in Red color



Replaces  
the 6301

## MDS3024B14+MD30-EMG1

One of the best medical emergency carts at a great crash cart price, this reliable, economical cart includes the MD30-EMG1 package with common accessories needed for most medical crash carts. Also includes Tente® casters, top quality drawer slides and a breakaway bar with 100 seals.



available in Red color

# EMERGENCY CRASH CARTS

Replaces the  
AL809B6-EMG

## MPA3027B06+MD30-EMG1

This lightweight aluminum cart comes with Harloff's MD30-EMG1 emergency accessory package which includes an oxygen tank holder, aluminum drop shelf, two utility hooks, two sets of standard 3" drawer dividers, a cardiac board with brackets and an adjustable IV pole. This cart comes with pre-threaded holes drilled into the cart so that each accessory can easily be installed with just a screwdriver upon delivery.

This lightweight break-away crash cart also features 5" ball-bearing Tente casters, a white antimicrobial plastic top, new full extension self-closing drawer slides, ergonomic side push handles and a fully-curved protective bumper. Help to ensure that your emergency care staff have all the supplies they need to save lives!



## MPA3030B06

Constructed from lightweight aluminum, this easy to maneuver crash cart uses self-closing, steel ball bearing, full extension drawer slides. Also features sleek, longer drawer pulls and premium 5" Tente® casters.



## MDS3024B05+MD30-EMG1

This short crash cart comes with the MD30-EMG1 emergency accessory package which includes an O2 holder, IV pole, cardiac board, drop shelf, utility hooks and drawer dividers. Comes with three 3" drawers, one 6" drawer and one 9" drawer.



Replaces the 6331

Replaces  
the 6416B

## MDS3030B06

Comes with a heat sealed powder coat steel breakaway lock bar with 100 seals. Also features 30" of vertical drawer space, an antimicrobial plastic top, protective bumpers, ergonomic side handles and more.



available in Red color



# EMERGENCY CRASH CARTS

MDS3030B16+MD30-FDR12

Similar to the 6403

## MDS3030B16

Replaces the 6400

Used in hospitals nationwide, this crash cart is ideal for emergency preparedness. Features include premium Tente® casters, full extension drawer slides, protective bumpers, an antimicrobial plastic top, ergonomic side handles and a breakaway locking bar with 100 seals.



## MDS3024B14

Replaces the 6300

This short emergency cart includes a heat sealed powder coat steel break-away lock bar w/100 numbered seals, 24" of vertical drawer space, full extension drawer slides, an antimicrobial and removable plastic top, protective bumpers, ergonomic handles and premium 5" quiet rolling Tente® casters.



This five-drawer cart has a bottom compartment with a 12" flip up door for ideal storage of large equipment. Other features include a side drop shelf, O2 holder, IV pole, utility hooks, premium casters, full extension drawer slides, an antimicrobial plastic top and a breakaway bar with 100 seals.



## MD30-EMG ACCESSORY PACKAGE

For MDS30/MPA30 carts, this Basic Package includes:

- Oxygen tank holder (O2HLDR)
- Cardiac board and brackets (MD-CARDBRD)
- IV pole (MD-IV2)



## MDS3030B06 +MD30-EMG

Provides an optimal crash cart setup for your facility with premium Tente® casters, top quality drawer slides and the MD30-EMG package which includes the most commonly requested accessories.



## MDS3024B05

Replaces the 6330

This M-Series crash cart features premium Tente® casters, full extension drawer slides and a breakaway bar with 100 plastic seals. M-Series carts allow numerous customizable options, including lock style, color, drawer configurations and even other heights and widths if you choose.



## M3DSI824B05+MD18-EMG

Narrow emergency carts like this can be situated in small, tight spaces while still storing an abundance of emergency supplies.

Comes with an MD18-EMG emergency accessory package which consists of an IV pole, cardiac board, side push handle and O2 holder.

Replaces the 3145B-EMG



## EMERGENCY CRASH CARTS

### MDS2430B06+MD24-EMG

This tall, medium width M-Series cart comes with Harloff's MD24-EMG emergency accessory package which includes an oxygen tank holder, cardiac board and IV pole.



### MD24-EMG ACCESSORY PACKAGE

**For MDS24/MPA24 carts, this Basic Package includes:**

- Oxygen tank holder (O2HLDR)
- Cardiac board and brackets (MD-CARDBRD)
- IV pole (MD-IV2)



### MD18-EMG ACCESSORY PACKAGE

**For MDS18/MPA18 carts, this Basic Package includes:**

- Oxygen tank holder (O2HLDR)
- Cardiac board and brackets (MD-CARDBRD)
- IV pole (MD-IV2)
- Aluminum Push Side Rail (SIDERAIL16)



### M3DSI824B05

Replaces the 3145B

Narrow emergency carts like this can be situated in small, tight spaces while still storing tons of supplies. This cart offers 24" of vertical drawer space, a painted breakaway locking bar with 100 plastic seals, a white antimicrobial plastic top and 3" stem casters.



### MDSI830B06+MD18-EMG

This mini-width M-Series cart comes with the MD18-EMG emergency accessory package which consists of a medical-grade IV pole, cardiac board, side push handle and oxygen tank holder.

Replaces the 3156B-EMG



# VALUE EMERGENCY CRASH CARTS

## V30-6EMG2

V-Series crash carts feature Harloff's quality construction offered at a value price. The accessory package includes the most requested exterior mounted accessories.



## V24-4EMG2

This crash cart is designed for healthcare facilities under budget pressure that require a durable, fully equipped emergency cart offering quality construction at a value price. Includes an IV pole, cardiac board, drop shelf, O2 holder, painted steel top rail, plastic top and drawer trays.



## V30-6EMG

V-Series crash carts feature Harloff's quality construction offered at a value price. This model includes a cardiac board, side mounted drop shelf, IV pole and O2 holder. This six drawer crash cart comes equipped with a breakaway bar with 100 plastic seals.



## V24-4EMG

V-Series carts provide a durable, fully equipped emergency cart offering Harloff's quality construction at a value price. Includes two 3" drawers, one 6" drawer and one 12" drawer plus an emergency accessory package.



## V24-4B

Designed for healthcare facilities under budget pressure, V-Series carts provide a durable emergency cart that offers Harloff's quality construction at a value price. This four drawer cart comes with a side rail, protective bumpers and 5" ball bearing casters.



## V30-6B

V-Series carts are designed for healthcare facilities with limited budgets. This crash cart has a top rail, breakaway bar and six storage drawers.





## MDS3030B09PED-EMG

The included IV pole, oxygen tank holder, cardiac board, drop shelf and utility hooks ensure you have all the necessary emergency equipment to treat pediatric patients. This M-Series cart is also highly customizable.



Replaces the 6401PEC

## PEDIATRIC EMERGENCY CARTS

### MDS3030B09PED

Designed for pediatric emergencies, this cart uses a multi-colored drawer labeling system comparable to the Broselow color coding method.



Replaces the 6400PEC

### MPA3030B09PED

This lightweight aluminum version of our pediatric cart has individual breakaway bars on each of its nine drawers, making restocking easy.



## ADDITIONAL MODELS

All carts can be equipped with any lock type and most can accommodate additional accessory packages.

Item #	Description
MDS3027B05	M-Series Steel Crash Cart, Medium Height, Standard Width, 5-Drawer, Breakaway Lock
MDS3027B06	M-Series Steel Crash Cart, Medium Height, Standard Width, 6-Drawer, Breakaway Lock
MDS3021B14	M-Series Steel Crash Cart, X-Short Height, Standard Width, 4-Drawer, Breakaway Lock
MDS3021B15	M-Series Steel Crash Cart, X-Short Height, Standard Width, 5-Drawer, Breakaway Lock
MDS2424B05	M-Series Steel Crash Cart, Short Height, Medium Width, 5-Drawer, Breakaway Lock
MPA3030B07	A-Series Tall, Standard Width, 7-Drawer cart with Breakaway Lock replaces AL810B7
MPA3027K06	A-Series Medium Height, Standard Width, 6-Drawer cart with Key Lock replaces AL809K6
MPA3027B06	A-Series Medium Height, Standard Width, 6-Drawer, Breakaway Lock replaces AL809B6

## CRASH CARTS ORDER SHEET

Item #	Description
MPA3027B05	A-Series Medium Height, Standard Width, 5-Drawer, Breakaway Lock replaces AL809B5
MPA3024B05	A-Series Short, Standard Width, 5-Drawer cart with Breakaway Lock replaces AL808B5
MPA2424B05	A-Series Aluminum Crash Cart, Short Height, Medium Width, 5-Drawer, Breakaway Lock
MPA1830B06	A-Series Narrow, Short Height, 6-Drawer cart with Breakaway Lock replaces AL3256B
MPA1824B05	A-Series Narrow, Short Height, 5-Drawer cart with Breakaway Lock replaces AL3245B
MDS3030B09PED+MD30-EMG3 MDS3030B09PED pediatric cart plus a MD30-EMG3 emergency accessory package	
V24-3B	V-Series 3-Drawer Short Crash Cart
V24-5B	V-Series 5-Drawer Short Crash Cart
V30-5B	V-Series 5-Drawer Tall Crash Cart

# CUSTOM AND UNIQUE PROJECTS

*Harloff has an unparalleled ability to do custom and unique projects. Let us design the perfect cart or storage system for your facility! We have helped to pioneer small cart innovations such as bayonet mounted ball-bearing drawer slides, fully modular customizable designs and carts that ship fully assembled.*

*The carts on the right are custom solutions requested by customers with options like independent breakaway locking bars, locking bars with antimicrobial white paint and undercarriage step stool holders. Step stools allow shorter nurses to gain additional height for effective chest compressions. Harloff can make custom changes to your casters, bumpers, locks, colors, drawer options, door options, top, height, width, depth, accessories and many other cart features. Allow us to help you create your cart just the way you want it!*



## CART OPTIONS

*Harloff has standard sizes, but our manufacturing capabilities allow us to build custom carts to your specifications. This cart is a double height, double width solution designed for a multi-facility off-site SPD department.*



## SECURITY OPTIONS

*Harloff offers multiple locking options, including advanced electronic solutions. This cart uses audible anti-theft and door ajar alarms, wireless access tracking and independent drawer locking options.*



## DOOR OPTIONS

*Our medical storage cabinets come with your choice of solid, glass-paneled or tambour doors, but sometimes custom options are needed. This two column MedStor Max cabinet was built with dual sliding doors secured with a key lock.*







## TOP QUALITY LOCKS FOR MAXIMUM SECURITY AND STORAGE

*Harloff manufactures the highest quality anesthesia storage carts in a variety of sizes, ideal for use in surgical environments to increase operating room efficiency and improve the patient experience. Anesthesia carts require proper security and strict access control to protect medications and supplies stored within. Many anesthesiologists maximize their cart's storage capabilities by outfitting them with back rails, tilts bins, side rails, hooks and other accessories, all of which are offered by Harloff.*

## ANESTHESIA CARTS

*Harloff anesthesia carts provide secure and convenient access to controlled substances, medications and other supplies that practitioners require. We offer multiple configurations of anesthesia carts for a variety of medical settings, available in both all-steel construction or light-weight aluminum with heat sealed powder coat finishes. Best of all, our anesthesia carts all come with Harloff's twelve year manufacturer's warranty.*



Engineered and  
Manufactured in the U.S.

MPA3030WLP16



# ANESTHESIA CARTS

Equip your cart with an accessory package for an ideal anesthesia cart setup, or select from other options or accessories! M-Series and A-Series offer you choices.

## MDS3030E06+MD30-ANS3

Replaces  
the 7456E

This popular M-Series anesthesia cart includes the MD30-ANS3 Anesthesia Package, Harloff's most comprehensive and robust package for anesthesiologists. It is constructed with full extension drawer slides, premium casters, protective bumpers, ergonomic handles, an antimicrobial plastic top and 30" of vertical drawer space. This complete setup offers an ideal cart, accessory and drawer configuration but, like any M-Series cart, can be customized to match your exact specifications.



## MD30-ANS3 ACCESSORY PACKAGE

**For MDS30/MPA30 carts\*, this Deluxe Anesthesia Package includes:**

- 3-Tier Back Rail (MD30-RBR3)
- (3) 12" Side Rails (SIDERAIL12)
- Large Utility Basket (BASKET3RC)
- Tape and Label Dispenser (LABELDISP)
- Catheter Box (CATHBOX18RC)
- (2) Standard 3"



Drawer Dividers with Tray (MD30-TRAYDIV3-S)

- (1) Basic 6" Drawer Divider (MD30-DIV6-B)
- 5-Compartment Tilt Bin w/Rail Clips (TILTBIN5RC)
- Sharps Container w/Wire Holder (MD-SHARPS-RC)
- Plastic Waste Container (MD-3GALWASTEDM)

## MD30-ANS ACCESSORY PACKAGE

**For MDS30/MPA30 carts, this Basic Anesthesia Package includes:**

- 3-Tier Back Rail (MD30-RBR3)
- (3) 12" Side Rails (SIDERAIL12)
- 5-Compartment Tilt Bin with Rail Clips (TILTBIN5RC)



## MPA3030E07+MD30-ANS

A light seven drawer e-lock anesthesia cart comes with a basic accessory package.



Replaces the AL810E7-ANS

## MPA3030EKC06

Replaces the  
AL810EKC6

This A-Series lightweight aluminum cart comes with ergonomic handles, fully curved bumpers, 5" quiet rolling Tente® casters and an electronic keypad lock with automatic re-locking drawers and access rights for up to 250 users.



## MDS3030K16+MD30-AIRWAYPKG

Replaces  
the 6457

This M-Series difficult airway cart securely stores laryngoscope blades, forceps, tubes, catheters and supplies to help manage intubations or failed airways. The included accessory package (MD30-AIRWAYPKG) includes a four-capacity side scope cabinet, an adjustable overhead shelf with tilt bins, utility hooks, a side rail and an oxygen tank holder. This six drawer key locking cart features an antimicrobial plastic top, highly mobile premium Tente® casters, an easily cleaned heat sealed powder coat finish and fully customizable drawer configurations to match all your needs.



## ANESTHESIA CARTS

### MD30-AIRWAYPKG ACCESSORY PACKAGE

**For MDS30 & MPA30 carts, this Difficult Airway Package includes:**

- Adjustable Overhead Shelf (MD30-ADJSHLF2)
- 5-Compartment Tilt Bin w/Rail Clips (TILTBIN5RC)
- Oxygen tank holder (O2HLDR)
- 4-Scope Capacity Side Scope Cabinet (680350)
- (1) 16" Side Rail (SIDERAIL16)
- (2) Utility hooks (UHOOKS2)
- (2) Standard 3" Drawer Dividers w/Tray (MD30-TRAYDIV3-S)
- (1) Basic 6" Drawer Divider (MD30-DIV6-B)



### MPA3027K06+MD30-ANS

This lightweight A-Series six drawer cart features our MD30-ANS anesthesia accessory package.



Replaces the AL809K6-ANS

### MDS3030K06

Customizable M-Series cart features:

- Key lock with two keys
- Full extension drawer slides
- Antimicrobial plastic top, removable
- 30" of vertical drawer space
- 5" quiet rolling Tente® casters



available in Sand color



# ANESTHESIA CARTS

MDS3030K06+MD30-ADJSHLF2-XBAR

## MPA3030WLP16

This A-Series cart has a premium CompX® wifi-capable lock with alarms and reports via email, text, fax and voice. It allows up to 3,000 user codes with a wireless audit of the last 15,000 access attempts via the control pad or proximity reader.



## MDS3030E16

This six drawer M-Series cart is secured with an electronic pushbutton lock with key lock override and factory-set user/supervisor codes.



Replaces the 7450E

This key locking M-Series cart can hold a printer or other light equipment with two adjustable overhead shelves. Each shelf measures approximately 14" x 27" and can hold up to thirty pounds. *Replaces the 6650*



## MDS3030K06+MD30-ANS2

Store the tools you need for putting patients safely under anesthesia with this M-Series cart. Industry best 5" quiet rolling Tente® casters provide excellent mobility between rooms while ergonomic side handles and protective bumpers ensure easy transport. Full extension drawer slides give easy access to your equipment, and this optimal drawer configuration gives you plenty of storage space. This cart includes Harloff's MD30-ANS2 Standard Anesthesia Package.



## MD30-ANS2 ACCESSORY PACKAGE

**For MDS30/MPA30 carts\*, this Standard Anesthesia Package includes:**

- 3-Tier Back Rail (MD30-RBR3)
- Sharps Container w/Wire Holder (MD-SHARPS-RC)
- (3) 12" Side Rails (SIDERAIL12)
- (2) Utility hooks (UHOOKS2)
- Catheter Box w/Rail Clip (CATHBOX18RC)
- 5-Compartment Tilt Bin w/Rail Clips (TILTBIN5RC)
- (1) Standard 3" Drawer Divider with Tray (MD30-TRAYDIV3-S)
- (1) Basic 6" Drawer Divider (MD30-DIV6-B)





## MDS3030K06+MD30-ANS

Replaces  
the 6451

This tall M-Series cart includes the MD30-ANS accessory package which features our most commonly requested items for anesthesia carts: a 3-tier-back rail, tilt bins and side rails for attaching rail clips.



Replaces  
the AL810E7

## MPA3030E07

A-Series carts are constructed from light, durable aluminum that allows you to effortlessly push them from one room to another. This 7-drawer cart has 30" of vertical drawer space plus a pushbutton e-lock.



# ANESTHESIA CARTS

## MDS3036K08

This M-Series XX-Tall high capacity cart is our tallest model at nearly 49" and provides a full 36" of vertical drawer space. Built for optimal storage but still provides high mobility.



## MPA3027E06+MD30-ANS

This A-Series six drawer cart utilizes our MD30-ANS accessory package and a pushbutton e-lock.



## MDS2430E06+MD24-ANS

This cart has all the standard M-Series bells and whistles plus the MD24-ANS Anesthesia Package.



## MD24-ANS ACCESSORY PACKAGE

**For MDS24/MPA24 carts, this Basic Anesthesia Package includes:**

- 3-Tier Back Rail (MD24-RBR3)
- (3) 12" Side Rails (SIDERAIL12)
- 5-Compartment Tilt Bin with Rail Clips (TILTBIN5RC)



# ANESTHESIA CARTS

## MD18-ANS ACCESSORY PACKAGE

**For MDS18/MPA18 carts, this Basic Anesthesia Package includes:**

- 3-Tier Back Rail (MD18-RBR3)
- (3) 12" Side Rails (SIDERAIL12)
- 5-Compartment Tilt Bin with Rail Clips (TILTBIN5RC)



## MDS1830E06+MD18-ANS

*This M-Series cart fits into tight areas that don't have much free space. Comes with the MD18-ANS accessory package consisting of back rails, side rails and tilt bins.*

*Replaces the AL3256E-ANS*



## MPA1830ELP06

*The ELP lock of this narrow A-Series cart comes with a keypad and proximity reader that enables you to use existing credentials to provide HID compatible keycard access for up to 250 users per cart. This lock also allows you to program the duration before an automatic cart re-lock, and to add, edit or delete users.*



*Replaces the AL3256ELP*

## MPA1830E06

*This A-Series aluminum lightweight cart comes with self-closing, full extension drawer slides and a push-button e-lock.*



*Replaces the AL3256E*

## MDS1830E06

*Replaces the 3156E*

*Narrow enough to fit in tight spaces but still providing a full 30" of vertical drawer space, this M-Series cart is tall enough to hold all your supplies and comes with an electronic pushbutton lock and plenty of customization options.*



## MPA1824K05+MD18-ANS

*This A-Series cart includes the MD18-ANS accessory package with a five compartment tilt bin, back rail and three 12" side rails.*

*Replaces the AL3245K-ANS*



# HARLOFF CAN CREATE CUSTOM SOLUTIONS FOR YOU: A CASE STUDY

*Mercy Medical Center approached Harloff to build an anesthesia cart with top notch security. Our solution included the following features:*

- A LockView® 5 compatible wireless lock utilizing 802.11g wireless protocol or Ethernet communication on existing networks. Used common security encryption methods; can program up to 3,000 user codes and audit the last 15,000 access attempts
- Built-in anti-theft drawer alarm
- Drawer ajar alarm that sounds after 30 seconds of inactivity with any drawer left open
- Able to use existing credentials and add new users with existing HID Prox and iClass compatible badges

*This cart also included an overhead shelf, 30" of vertical drawer space and standard M-Series cart features. A-Series and M-Series carts are highly customizable and offer multiple locking solutions – let us build the perfect cart for you!*



## MDS1824K05

*Replaces  
the 3245K*

*This M-Series mini width cart has tons of storage space but has a small enough profile to fit in tight areas. This highly mobile cart comes with 5" quiet rolling Tente® casters, stabilizing pontoon bumpers and smooth-gliding full extension drawer slides.*



# ANESTHESIA CARTS

## MPA1830K06+MD18-ANS

*Includes Harloff's MD18-ANS accessory package consisting of a back rail, side rails and a 5-compartment tilt bin. This A-Series cart's light aluminum build and narrow cabinet fit well into tight spaces.*



*Replaces the  
AL3256K-ANS*

## MPA1830EKC06

*This A-Series cart has an auto-locking electronic keypad that grants access rights for up to 250 users. It also provides 30" of vertical drawer space, full extension self-closing drawer slides and premium Tente® casters.*



*Replaces the  
AL3256EKC*

## MPA1824E05+MD18-ANS

*This lightweight narrow anesthesia cart comes with a three tiered back rail, side rails and a five compartment tilt bin as well as self-closing drawer slides.*



*Replaces the  
AL3245E-ANS*



# ANESTHESIA CARTS

M3DS1830K06+MD18-ANS

## M3DS1824E05

Replaces  
the 3145E

This narrow M-Series cart packs a lot of storage space into a small frame. It features full-extension drawer slides, an antimicrobial white plastic top, 3" stem casters, 24" of vertical drawer space and customizable drawer configurations. Cart contents are kept well-secured by an electronic pushbutton lock that accepts up to 20 different user codes.



## MPA1824E05

Replaces  
the AL3245E

This narrow aluminum A-Series cart uses a pushbutton lock which provides code access for up to 19 different users and one supervisor. The narrow profile of this 5-drawer cart gives it a minimal footprint, perfect for facilities where every inch of floor space is at a premium.



Includes Harloff's MD18-ANS basic anesthesia accessory package and 3" stem casters.

Other features include full extension drawer slides, an antimicrobial plastic top and a key lock.



## ADDITIONAL MODELS

Item #	Description
All carts can be equipped with any lock type and most can accommodate additional accessory packages.	
MDS3027K05	Medium 5-Drawer M-Series + Key Lock
MDS2430K15	Tall 5-Drawer Medium Width M-Series
MDS1830EKC06	Tall M-Series Mini + Elec. Keypad Lock
MDS1824EKC05	Short M-Series Mini + Elec. Keypad Lock
MPA3030K06	This Tall 6-Drawer A-Series is a popular all-around model used w/accessory packages
MPA3030K16	An ideal model for difficult airway carts
MPA3030E06	Tall A-Series replaces the AL810E6
MPA3030E07	Tall A-Series replaces the AL810E7
MPA3027E06	Medium 6-Drawer cart replaces AL808E5
MPA3024E05	Short 5-Drawer A-Series + E-Lock

Customize any A-series or M-series cart with your choice of accessory packages to create the ideal anesthesiologist room setup.



Item #	Description
MPA2427E05	Medium Height & Width 5-Drawer A-Series
MPA1830K06	Tall 6-Drawer A-Series replaces AL3256K
MPA1830E06	Tall A-Series Mini replaces the AL3256E
MPA1824K05	Short 5-Drawer A-Series replaces AL3245K
MPA1824EKC05	Short A-Series Mini replaces AL3245EKC
MPA1821ELP03	X-Short 3-Drawer A-Series Mini + Elec. Keypad Lock with Proximity Reader replaces the AL3233ELP

# IT TOOK A PANDEMIC TO REALIZE THE VALUE OF INFECTION CONTROL

*The need for proper infection control and isolation measures at medical facilities across the globe has never been higher than today. Having a secure place to store PPE for your hard-working staff is vital. Harloff isolation carts typically feature deep drawers and removable antimicrobial tops. Multiple sizes, configurations and colors of carts are available although most come in a yellow color indicative of isolation areas.*

## ISOLATION CARTS

***Isolation carts, also referred to as infection control carts, are intended to help prevent the spread of infectious illness and disease within a healthcare setting. Storage solutions for PPE supplies help to reduce the risk of pathogen transmission between patients, healthcare workers and visitors.***

***Stopping the spread of VREs, MRSAs and other infectious microbes can help save hospitals millions of dollars as well as improve patient outcomes. Our carts help lower hospital costs and liability by limiting the number of healthcare acquired infections.***



Engineered and  
Manufactured in the U.S.



MPA2430EKC16

# INFECTION CONTROL CARTS

## MDS2424E05

Replaces  
the 4245E

Formerly called the Mini24, this cart with 18.5" wide drawers is a happy medium between the full width cart and the minis. Includes full extension drawer slides, fully curved protective bumpers, an anti-microbial plastic top and an electronic pushbutton lock with key lock override.



## M3DSI830K04-PPE

This four drawer, key locking PPE cart has ample storage for face shields, masks, medical gloves, gowns and footwear plus it comes with a Sani-Wipe holder and triple glove box holder. It uses one of our most popular cart bodies.



**QuickShip**

available in Yellow color

## MDS3030K03

Replaces  
the 6510

With one 6" drawer and two 12" deep drawers, this cart has plenty of storage space for all your PPE. Ergonomic side push handles and premium 5" Tente® casters provide optimal mobility and full extension drawer slides allow easy access to drawer contents.



## MDS3030K04

Replaces  
the 6520

This top selling infection control cart has one 3" drawer and three 9" deep drawers. Further customize your cart with accessories such as drawer dividers, tilt bins and glove boxes.



**QuickShip**

available in Yellow color

## M3DSI824K03

Replaces  
the 3143K

This short mini-width cart helps organize all your PPE and other supplies in one mobile, secure location. Three large drawers, each with full extension drawer slides offer up to 24" of vertical drawer space. Includes an antimicrobial plastic top, 3" stem casters and multiple customizable options.



## M3DSI830K04-VC

This vaccine cart comes with four drawers of dry storage for vaccine preparation, along with a triple glove box holder and sharps container for syringes. This is a great solution for vaccination efforts that don't require cold storage.





## VR4300K-AC

Administer the COVID-19 vaccine and other refrigerated vaccines easily with this vaccine cart, which has a built-in medical-grade lab quality refrigerator that keeps supplies at 2° C with 2.4 cubic feet of space for cold storage. Three top drawers provide space for additional supplies.



## INFECTION CONTROL CARTS

Replaces  
the AL3244K

### MPA1824K04

This narrow lightweight aluminum cart offers 24" of vertical drawer space, secured by a standard key lock with two keys. Features premium 5" Tente® casters, a white antimicrobial plastic top, new full-extension self-closing drawer slides and the industry's only fully-curved protective bumper.



Replaces  
the 3154K

### M3DS1830K04

One of our top-selling items during the pandemic, this cart combines premium PPE storage abilities with a small footprint in one economical package. Comes with a key lock, 3" single wheel stem casters and 30" of vertical drawer space.



### MPA2430EKC16

This lightweight aluminum cart has a highly secure electronic keypad lock plus a six drawer configuration for maximum storage. All A-Series carts include 5" premium casters, pontoon bumpers and self-closing drawer slides.



### 3123

The IV Start Cart is designed to reduce the risk of accidental needle sticks to patients and staff while improving the efficiency of IV and epidural procedures. Features a dual work surface, articulating arm, waste container and Sharps container. Comes with either a latch or key lock to secure the four full extension drawers.



### 3256KP8

This phlebotomy cart has 30" of vertical drawer space secured with either a key or latch lock and it comes with two 6-compartment tilt bins, a raised back rail, a waste container, a single glove box holder and extended push handle. Phlebotomy carts keep your supplies nearby, are mobile, provide a dedicated work surface for procedures and allow quick disposal of used needles.



# INFECTION CONTROL CARTS

## V30-4K

V-Series Carts are value-priced and used in a variety of hospital departments. This drawer configuration is a popular choice for isolation carts and features a key lock, side rail push handle and plenty of room for accessories.



## V30-3K

This tall three drawer V-Series cart comes with a key lock, side rail push handle, corner bumpers, 5" ball bearing casters and 30" of vertical drawer space with plenty of room for PPE and other supply storage.



## MDS3030E03

Replaces the 7510E

Add additional security to your cart with an electronic pushbutton lock that stores access codes for up to twenty users.

Other features include:

- Full extension drawer slides
- Antimicrobial plastic top, removable
- Ergonomic side push handles
- 30" of vertical drawer space
- Fully curved protective bumper
- 5" quiet rolling Tente® casters



## OTHER INFECTION CONTROL CARTS

Item #	Description
M3DS1824E03	M-Series Short Mini Cart, 3-Drawer, 3" Stem Casters, E-Lock, replaces 3143E
M3DS1824E04	M-Series Short Mini Cart, 4-Drawer, 3" Stem Casters, E-Lock, replaces 3144E
M3DS1824K04	M-Series Short Mini Cart, 4-Drawer, 3" Stem Casters, Key Lock, replaces 3144K
M3DS1830K16	M-Series Tall Lightweight Aluminum Cart, Mini, 6-Drawer, 3" Stem Casters, Key Lock
MDS2424E03	M-Series Short Medium Width Cart, 3-Drawer, E-Lock, replaces 4243E
V30-3K	V-Series 3-Drawer Cart, Tall, Key Lock
V24-3K	V-Series 3-Drawer Cart, Short, Key Lock
V30-4B	V-Series 4-Drawer Cart, Tall, Breakaway Lock
V30-4K	V-Series 4-Drawer Cart, Tall, Key Lock

## ADDITIONAL MODELS

### Vaccine Carts

VR4302K-AC Two drawer version of the VR4300K-AC

VF4250K-AC Two drawer vaccine freezer cart with a 3.2 cubic feet Accucold® freezer that cools down to -25° C

### Phlebotomy Carts

3254KPB Four drawer version of the 3256KPB

3256LPB 3256KPB with a latch lock

**INFECTION CONTROL CARTS  
HELP LOWER MEDICAL COSTS  
BY REDUCING THE NUMBER OF  
HOSPITAL ACQUIRED INFECTIONS**



# IMPROVE EXAM EFFICIENCY WITH TREATMENT & PROCEDURE CARTS

*Treatment carts are the perfect solution when you need to store topical ointments, eye, ear, nose and throat supplies, wound care supplies and other bulkier nursing supplies including surgical apparel, masks, dressing disposal bags and blankets. These carts are seen throughout medical facilities including in hospital floors, nursing units, specialty clinics, operating rooms and treatment rooms.*

M3DS1830K06

## CARTS FOR TREATMENT/ PROCEDURES

*Treatment and procedure carts are the perfect solution for storing all the supplies for wound care and other general patient needs in one place. Our carts help keep numerous small supplies well organized for efficient and effective application to improve patient outcomes and experience.*

*Harloff offers a variety of specialized cart configurations to fit your needs, including top of the line A-Series lightweight aluminum carts, M-Series all-steel construction carts and economy lines like ETC carts or V-Series models.*



Engineered and  
Manufactured in the U.S.





# TREATMENT/PROCEDURE CARTS

## M3DS1830K06

Replaces  
the 3156K

This narrow procedure cart can fit in tight spaces and features elegant curved drawer pulls, full-extension drawer slides, 30" of vertical drawer space and 3" stem casters.



## MDS3030K02DR

This key locking multi treatment cart includes one 3" drawer, one 6" drawer and one 21" storage compartment with an adjustable shelf secured behind a double door with a separate lock.



## MPA3030E06

This lightweight A-Series aluminum cart is equipped with industry best 5" quiet rolling Tente® casters, ergonomic side push handles and fully curved protective bumpers. Also features curved drawer pulls for a durable and elegant look, self-closing full extension drawer slides and 30" of customizable vertical drawer space. This cart comes equipped with an electronic pushbutton lock.



## MDS3030E16

Replaces  
the 7450E

This six drawer cart has full extension, smooth gliding drawer slides, a push-button e-lock with key lock override, ergonomic push handles, industry best Tente® casters and your choice of drawer configurations and color.



## MDS3024E14

Replaces  
the 7350E

This short four-drawer cart comes with a keyless electronic pushbutton lock that stores up to 20 user codes and it can cycle up to 75,000 times before power depletion with standard batteries. It also comes with industry best Tente® casters, full extension drawer slides and ergonomic push handles.



## M3DS1824K05

This narrow key locking procedure cart has a small footprint making it ideal for overcrowded rooms with limited space.

Replaces the 3145K



Replaces  
the 6450

## MDS3030K16

Built to last, this cart has a sturdy, uni-body steel construction with a scratch resistant heat sealed powder coat finish, inside and out, in your choice of over 20 colors. It also features premium Tente® casters, the industry's only fully curved bumper, full extension drawer slides and a removable antimicrobial top. Enhance your cart with one of our accessory options.



Replaces  
the AL808E5

## MPA3024E05

This A-Series cart features a light-weight aluminum construction, new curved drawer pulls, self-closing drawer slides, protective bumpers and much more. It comes configured with an electronic pushbutton lock, but users can choose between 10 different locking options.



## V30-6K

V-Series carts are value priced for facilities with limited budgets. This tall six-drawer cart features a key lock with two keys, an aluminum side rail/push handle and 30" of vertical drawer space.



## ETC-4

This four drawer ETC cart utilizes a key lock with two keys, a removable plastic top, drawer fronts with integrated pull handles, roller drawer slides and 3" stem casters. ETC carts are budget friendly, basic carts ideal for bedside use.



## V24-5K

This short five drawer V-Series cart has 24" of vertical storage using three 3" drawers, one 6" drawer and one 9" drawer. Comes with a key lock, side rail/push handle, corner bumpers, 5" ball-bearing casters and full extension drawer slides.



## ETC-5

ETC carts are an ideal solution for budget driven facilities requiring mobile storage for supplies. The short height of these carts is perfect for bedside use, but their compact size can fit clinics, doctor offices, schools and veterinary centers too. This key locking cart features 3" stem casters.





# TREATMENT/PROCEDURE CARTS

## MDS3024K14

Replaces  
the 6350

This short four drawer M-Series cart features industry best Tente® casters for quick movement, full extension drawer slides and ergonomic handles. All M-Series carts allow numerous customizable options including lock style, color, drawer configurations and even custom heights and widths if needed.



## MPA2427E05

This e-locking A-Series cart is made from extremely lightweight aluminum, ideal for carts that are moved frequently. Features curved drawer pulls and self-closing full extension drawer slides.



## V30-2K



## V24-3K

## ADDITIONAL MODELS

Item #	Description
MPA3030K06	Tall A-Series with 3,3,3,6,6,9 drawer config
MPA3030K07	Tall A-Series with 3,3,3,3,3,6,9 drawer config
MPA3024K05	Short A-Series replaces the AL808K5
MPA1830E06	Tall A-Series Mini replaces the AL3256E
MPA1824EKC05	Short A-Series replaces the AL3245EKC
MPA1821ELP03	X-Short A-Series replaces the AL3233ELP
MPA1821K03	X-Short A-Series replaces the AL3233K
MDS3030K06	MDS3030K16 with 3,3,3,6,6,9 drawer config
MDS3030EKC06	MDS3030K06 + Electronic Keypad Lock
MDS3036K08	XX-Tall M-Series with 36" of vertical drawers
MDS1821E14	X-Short M-Series with four drawers + E-lock
V24-4K	V-Series 4-Drawer Short Cart with Key Lock
V24-5K	V-Series 5-Drawer Short Cart with Key Lock

**ETC AND V-SERIES VALUE LINE CARTS ARE COST EFFECTIVE AND OUR SPACE SAVER DESIGNS CAN BE USED IN AREAS WITH LIMITED SPACE TO SECURE SUPPLIES**

Item #	Description
V24-8K	V-Series 8-Drawer Short Cart with Key Lock
V30-2K	V-Series 2-Drawer Tall Cart with Key Lock and 21" Locking Storage Compartment with Shelf
V30-5K	V-Series 5-Drawer Tall Cart with Key Lock
V30-6K	V-Series 6-Drawer Tall Cart with Key Lock
ETC-3	Economy 3-Drawer Short Cart with Key Lock, plus a 12" Bottom Open Storage Compartment
ETC-1-1D	Economy 1-Drawer Short Cart with Key Lock, plus a 17" Bottom Open Storage Compartment with Shelf